LUX-SUTURES 22 Gruuss-Strooss, 9991 Weiswampach, Luxemburg	DATA SHEET PROTECTIVE MASKS FFP2	F-Q-143E/R01
		2020 / 1
		Page 1 / 12

Content:

1. Data Sheet

- a) Description
- b) Filtration capacity
- c) GRAPHENE – Performance compared to other masks
- d) Three main advantages
- e) Article number - Packaging
- f) Practical use
- g) Use limitation
- h) Information on storage

2. Test reports

1. Datasheet

- a) Description

LUX-SUTURES 22 Gruuss-Strooss, 9991 Weiswampach, Luxemburg	DATA SHEET PROTECTIVE MASKS FFP2	F-Q-143E/R01
		2020 / 1
		Page 2 / 12

- **FFP2 - Graphene - protective mask.**
- **This is a new type of face mask whose main filter layer consists of a polypropylene graphene material which is applied to the melt blown fabric. This is the main filter layer of the mask.**
- **Graphene gives the masks stronger antibacterial properties, better air permeability and longer application time. The graphene mask uses the nano blade effect of the graphene material to destroy the cell wall of bacteria.**
- **Compared to the melt blown material used in conventional masks, the new graphene polypropylene layer has a unidirectional permeability function, making the mask more comfortable for users.**
- **The two ear loops can be adjusted in length with a small clip to ensure a perfect fit.**
- **A foam covers on the noseband improve the wearing comfort.**

b) Filtration capacity

LUX-SUTURES 22 Gruuss-Strooss, 9991 Weiswampach, Luxemburg	DATA SHEET PROTECTIVE MASKS FFP2	F-Q-143E/R01
		2020 / 1
		Page 3 / 12

- The filtering system of this mask is 99% effective against particles larger than 0.3 microns and 80% effective against anything smaller than 0.3 microns.
- The graphene is a "laser-induced graphene", a microporous conductive foam that can fix bacteria.
- Unlike other FFP2-KN95 masks, which fix bacteria and viruses to the outside of the mask, graphene masks can kill these particles.
- The size of a bacterium is 300nm to 1000nm and that of a virus 100-150nm. According to tests, the graphene mask is able to retain these small viruses.
- The BFE (Bacterial Filtering Efficiency) is > than 99%.
- The VFE (Virus Filtering Efficiency) is > than 99.97% - 99.98%.

c) GRAPHENE – Performance compared to other masks

- The graphene face mask makes use of the nano-knife effect of the graphene material to destroy the cell wall of bacteria.
- Graphene is a two-dimensional crystalline material, an allotropic form of carbon whose stack forms graphite, a structure 200 times harder than steel.
- Thanks to innovative technology, it becomes three-dimensional, which corresponds to a sharp, rounded blade with a thickness of 1 to 2 nm, capable of killing bacteria and viruses.
- The effectiveness of a graphene mask is 30 times greater than that of a normal face mask and 10 times greater than that of KN95

Graphene, as a two-dimensional carbon nanomaterial, can form a natural 'nano knife' on melt-blown cloth, which, under the action of natural respiration, destroys the cell walls of germs and plays a good bactericidal role.

d) Three main advantages:

- 1. Much higher efficiency: Nano-high-tech material is better than the standard N95 or FFP2**
- 2. Longer use - up to 200 hours while maintaining form and function with 99% BFE.**
- 3. Much more comfortable - easier to wear due to better fit, good air permeability, does not fog up when wearing glasses.**

e) Packaging:


**2 masks in a plastic bag
40 Masks / Box**

f) Practical use:

1. Wash hands or use clean gloves.
2. Unfold the mask completely.
3. Put on the mask.
4. Position ear loops correctly
5. Optimise the fit of the mask
6. Adjust the nose bridge: Press the nose clip lightly on the bridge of the nose with your index finger. With the other hand you adjust the nose clip to the shape of your nose. This is to ensure that there is no gap above the bridge of the nose and that the mask fits tightly.

To check that the mask fits correctly, place both hands on the mask and breathe deeply.

- If an airflow is felt in the nose area, readjust and tighten the nose clip.
- If you feel an airflow around the edges of the mask, reposition the mask and ear loop for one.

g) Use limitation

LUX-SUTURES 22 Gruuss-Strooss, 9991 Weiswampach, Luxemburg	DATA SHEET PROTECTIVE MASKS FFP2	F-Q-143E/R01
		2020 / 1
		Page 6 / 12

① Do not use the mask or enter or stay in a contaminated area under the following circumstance:

- a) Atmosphere contains less than 19.5% oxygen.
- b) If you smell or taste contaminant.
- c) For protection against gases or vapors.
- d) Contaminants or their concentrations are unknown or immediately dangerous to life or health.

② Do not modify or misuse the mask.

③ Do not use the mask with facial hair or any other conditions that may prevent a good face-seal, the requirements for leakage will be achieved.

④ Particle filtering half mask need to be inspected prior to each use to assure there are no holes in the breathing zone other than punctures around and staples and no damaged has occurred. Enlarged holes resulting from ripped or torn filter material around staple punctures are considered damage.

⑤ This mask helps protect against certain particulate contaminants but does not eliminate exposure to the risk or contracting disease or infection. Misuse may result in sickness or death.

h) Information on storage

LUX-SUTURES 22 Gruuss-Strooss, 9991 Weiswampach, Luxemburg	DATA SHEET PROTECTIVE MASKS FFP2	F-Q-143E/R01
		2020 / 1
		Page 7 / 12

- **Store the packaging away from light and moisture.**
- **Shelf life: 3 years from date of production.**

 <p>Storage Temperature Range: -20°C to + 30°C</p>	 <p>Storage Maximum Relative Humidity: <80%</p>	 <p>Shelf Life: 3 years from date of manufacture</p>	 <p>See information supplied by the manufacturer</p>
---	---	---	---

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